Efficient Brix Control for Beverages Simple, Rugged, and Trustworthy



Smart design With a vertical design and a sampling funnel, the Easy Bev refractometer simplifies workflows. The throughput is maximized while minimizing sedimentation on the prism when measuring samples such as juices.

Intuitive operation

Instrument set up is simplified, allowing first measurements within minutes: the intuitive app-based user interface minimizes user training and the integrated user guidance prevents operator errors.

Immediate QC results

Routine measurements are simplified and acceptable samples are instantly identified by setting up a range of tolerable limits. The color-coded results allow the user to immediately recognize if a measurement is within limits.

Easy performance verification

SimpleCheck[™] sucrose solutions for Brix performance verification ensure reliable results. The SimpleCheck[™] embedded workflow guides operators through the instruments' performance verification.



Easy Bev Refractometer

The Easy Bev Refractometer offers everything needed for daily routine applications: fast measurements, intuitive app based user interface, minimal operator dependency and more:

- Fast and reliable temperature control (can be turned off for instant results)
- LongClick[™] press and hold the home screen to start routine analyses in 1 second
- SimpleCheck[™] certified standards
- EasyDirect PC software for efficient data management





Technical and Functionalities Specifications

| EasyPlus Refractometers | | Easy Bev | |
|--|---|---|--|
| Refractive Index (nD) | Measuring range | 1.32–1.58 | |
| | Accuracy | ±0.0001 | |
| | Repeatability | ±0.00005 | |
| | Resolution | 0.0001 | |
| Brix (% w/w) | Measuring range | 0-100 | |
| | Accuracy | 0.051) | |
| | Repeatability | 0.05 | |
| | Resolution | 0.01 | |
| Temperature (°C) | Temperature range | 15–25 (59–77 °F) | |
| | Resolution | 0.01 | |
| | Temperature control | Yes ²) | |
| Measurement scales | Integrated concentration tables and functions | Refractive index, Brix, Brix compensated 20 °C, Inverted sugar, HFCS 42/55, Oechsle | |
| Number of methods | | 6 | |
| Internal results storage | | 203) | |
| Result export | | To USB stick (.csv) | |
| Color-coded results identification | | Yes | |
| SimpleCheck™ standards | | Yes | |
| Barcode reader connectivity | | Yes | |
| On-screen help text | | Yes | |
| Languages | | English, Chinese, German, Spanish, French, Italian, Portuguese, Russian, Bahasa Indonesia, Korean, Polish, Thai, Turkish | |
| Display | | 4.3" Capacitive touch screen | |
| Connectivity | | 3 × USB-A (keyboard, mouse, barcode reader, printer) | |
| | | 1 × USB-B (EasyDirect PC Software) | |
| PC Software | | EasyDirect Density and Refractometry PC Software | |
| | | (installation on a single PC, connect up to 3 density meters and/or refractometers) | |
| Minimal sample volume | | Approx. 50 mL | |
| Instrument dimensions (W × H × D) (mm) | | 198 × 382 × 200 | |
| | | $(7.7 \times 15 \times 7.8 \text{ inches})$ | |
| Service | | Professional Installation (StarterPac) | |
| | | Professional calibration | |
| Standard compliance | | www.mt.com/dere-norms | |
| ¹⁾ Accuracy of 0.07 at approx. 10% Brix and increasing up to 0.04 at approx. 80% Brix | | ³⁾ Unlimited results storage when transferred directly to a USB stick or to EasyDirect PC software | |

²⁾ Temperature control can be turned off for instant results

Instruments and Accessories

| Material Description | Material Number |
|--|-----------------|
| EasyPlus Refractometer – Easy Bev | 30468674 |
| EasyDirect Density and Refractometry PC Software | 30451628 |
| USB transfer cable EasyDirect PC Software | 51191926 |
| Printer USB-P25 | 11124301 |
| Barcode reader | 21901297 |
| USB cable for barcode reader | 21901309 |
| USB keyboard | 51192111 |
| SimpleCheck™ nD Water 1.333 @ 20 °C | 30467245 |

| Material Description | Material Number |
|--|-----------------|
| SimpleCheck™ nD Bromonaphthalene 1.657 @ 20 °C | 30467248 |
| SimpleCheck™ nD Dichlorotoluene 1.546 @ 20 °C | 30467247 |
| SimpleCheck™ nD Dodecane 1.421 @ 20 °C | 30467246 |
| SimpleCheck™ Brix Sucrose/water 10% @ 20 °C | 30467238 |
| SimpleCheck™ Brix Sucrose/water 12% @ 20 °C | 30467239 |
| SimpleCheck™ Brix Sucrose/water 20% @ 20 °C | 30467240 |
| Wall mount kit | 30465247 |
| Tubing set Easy Bev | 30472921 |

METTLER TOLEDO Group

Analytical Division Local contact: www.mt.com/contacts

www.mt.com/EasyPlus-Refractometry _

Subject to technical changes © 05/2018 METTLER TOLEDO. All rights reserved 30475527 Group MarCom 2470 RT/KA For more information